

SEQUENCE LISTING

<110> Dolnick, Bruce

<120> Compositions and Method for Identifying Thymidylate Synthase Modulators

<130> 03551.0143

<150> US 60/428,345

<151> 2002-11-22

<160> 9

<210> 1

<211> 150

<212> DNA

<213> artificial sequence

<220>

<223> TBE cassette

<400> 1

gataagcttc ctccgtcccc cgcccgccgc gccatgcctg tggccggctc
50
gtcagtcagg ctagctatag cggacttggg ccagtttat ggcttccagt
100
ggaggcattt tggggcagaa tacagagata tggaatcaga ttaagcttgc
150

<210> 2

<211> 42

<212> DNA

<213> artificial sequence

<220>

<223> TBE1 in the forward orientation

<400> 2

cctccgtccc ccgcccgcgc cgccatgcct gtggccggct cg
42

<210> 3

<211> 71

<212> DNA

<213> artificial sequence

<220>

<223> TBE2 in the forward orientation

<400> 3

ggacttgggc ccagtttatg gcttccagtg gaggcatttt ggggcagaat
50
acagagatat ggaatcagat t
71

<210> 4

<211> 42

<212> DNA

<213> artificial sequence

<220>

<223> TBE1 in the reverse orientation

<400> 4

gctcgcccg tgtccgtacc gcgcccgcg cccctgcct cc
42

<210> 5

<211> 71

<212> DNA

<213> artificial sequence

<220>

<223> TBE2 in the reverse orientation

<400> 5

ttagactaag gtagagagac ataagacggg gttttacgga ggtgaccttc
50
ggtatttgac ccgggttcag g
71

<210> 6

<211> 21

<212> DNA

<213> artificial sequence

<220>

<223> PCR sense primer for TS mRNA

<400> 6

tttggacagc ctgggattct c
21

<210> 7

<211> 21

<212> DNA

<213> artificial sequence

<220>

<223> PCR antisense primer for TS mRNA

<400> 7

aaagcaccct aaacagccat t

21

<210> 8

<211> 21

<212> DNA

<213> artificial sequence

<220>

<223> PCR sense primer for beta-actin

<400> 8

cagctcacca tggatgatga t

21

<210> 9

<211> 18

<212> DNA

<213> artificial sequence

<220>

<223> PCR antisense primer for beta-actin

<400> 9

ccagacgcag gatggcat

18

BFLODOCS: 873193 v1 (\$PRD01!.DOC)